

SURGICAL



## EQUIPMENT DRAPES

MICROSCOPE DRAPES

**ECOLAB**<sup>®</sup>

[www.ecolab.com](http://www.ecolab.com)

# MICROSCOPE DRAPES

Optical quality.  
Ease of use.  
Wide range.



# COVERING VIRTUALLY EVERY MICROSCOPE FOR MORE THAN TWO DECADES

Ecolab Surgical has worked in partnership with customers and leading manufactures to focus on providing you with superior imaging technology. Ecolab Surgical has developed and tested drapes for many years, and covers thousands of microscopes every day based on the premise that you need the best to do the best.

## MICROTEK Angled Glass Lens: One of the best in spectral reflectance (glare)

This test uses a spectrophotometer to measure the consistency of the light reflected off of the surface of a lens. If various wavelengths pass through the lens material differently, the image will be out of focus. A lower test result means a lens that will allow more light to pass through it.

Vari-Lens 2 and our patented Armatec™ range offers just that. Make sure you contact us for a free trial or a consult with one of our sales advisors.

- Very high scratch resistance
- Exceptional light transmission
- Reduced glare
- Enhanced definition (resolution)

Our Offer	Your Benefit
Optical excellence	Exceptional quality of vision
Wide range of covers	One high quality supplier for all your needs
Easy draping solution for virtually every microscope	First time right, quick draping solution, saving valuable time



Scan code for Microscope drape video

## OUR LENSES... AN EVOLUTION TO OPTICAL EXCELLENCE



### FLAT LENS

Polystyrene Lens providing sterile barrier and to cover basic optical and needs.



### VARI-LENS2

Angled Coated Glass Lens combining superior optical quality and cost effectiveness.



### CLEAR LENS

Curved Polycarbonate Lens providing sterile barrier and covering higher optical needs.



### ARMATEC

Innovative housing concept specially designed for Leica microscopes featuring a 'click' closure for secure attachment on the objective with an angled coated glass lens.



### VARI-LENS

Angled Coated Glass Lens featuring unique rotational option providing superior optical quality.

## RECOMMENDED MICROSCOPE DRAPES FOR ZEISS™ MD MICROSCOPES (65MM LENS)

Zeiss

Zeiss Microscope Body		Suggested Drapes		
Ophthalmology	Objective Lens Diameter (mm)	Flat Lens	ClearLens	Vari-Lens2
VISU	65	4841UK	4841CLEU	6100VL2
OPMI Lumera	65	4841UK	4841CLEU	6100VL2
OPMI Pico	65	4841UK	4841CLEU	6100VL2
ENT	Objective Lens Diameter (mm)	Flat Lens	ClearLens	Vari-Lens2
PRO ERGO	65	4842UK	4842CLEU	6100VL2
PRO MAGIS	65	4842UK	4842CLEU	6100VL2
OPMI Movena	65	4842UK	4842CLEU	6100VL2
OPMI Pico	65	4841UK	4841CLEU	6100VL2
OPMI Sensera	65	4862UK	4862UKCL	6100VL2
OPMI VARIO	65	4840UK	4840UKCL	6100VL2
OPMI 1 FC	65	4841UK	4841CLEU	6100VL2
TIVATO 700	65	4840UK	4840UKCL	6100VL2
Neurosurgery	Objective Lens Diameter (mm)	Flat Lens	ClearLens	Vari-Lens2
OPMI PENTERO	65	4842UK	4842CLEU	6120VL2
OPMI VARIO	65	4840UK	4840UKCL	6120VL2
TIVATO 700	65	4842UK	4842UKCL	6120VL2
Plastics	Objective Lens Diameter (mm)	Flat Lens	ClearLens	Vari-Lens2
OPMI PENTERO	65	4842UK	4842CLEU	6120VL2
OPMI VARIO	65	4840UK	4840UKCL	6100VL2















For more information about the drape sizes and packed quantities see page 5.

“I cannot afford to miss any details while using my microscope”




## ZEISS DRAPE SIZES AND PACKED QUANTITIES (65MM LENS)

Zeiss

Part. no	Sketch	Description	Qty/cs	Size
4861EU		Microscope drape for Zeiss™ MD with handle covers	20	76x163 cm
4831UK 4831UKCL		Microscope drape for Zeiss™ MD Microscope drape for Zeiss™ MD, featuring Clearlens™	20	102x61x163 cm
4841UK 4841CLEU		Microscope drape for Zeiss™ MD Microscope drape for Zeiss™ MD, featuring Clearlens™	20	117x163 cm
4862UK 4862UKCL		Microscope drape for Zeiss™ MD Microscope drape for Zeiss™ MD, featuring Clearlens™	10	117x203 cm
4840UK 4840UKCL		Microscope drape for Zeiss™ MD Microscope drape for Zeiss™ MD, featuring Clearlens™	10	117x267 cm
4842UK 4842CLEU		Microscope drape for Zeiss™ MD Microscope drape for Zeiss™ MD, featuring Clearlens™	10	137x267 cm
4838UK 4838CLEU		Microscope drape for Zeiss™ MD Microscope drape for Zeiss™ MD, featuring Clearlens™	10	117x305 cm
4865CLEU		Microscope drape for Zeiss™ MD, featuring Clearlens™	10	137x381 cm
6100EU		Vari-Lens™ microscope drape for Zeiss™, featuring rotational angled lens	10	122x305 cm
6100VL2		Microscope Drape, Zeiss™ with Vari-Lens2	10	122 x 305 cm
6120EU		Vari-Lens™ microscope drape for Zeiss™, featuring rotational angled lens	10	132 x 381 cm
6120VL2		Microscope Drape, Zeiss™ with Vari-Lens2	10	132 x 381 cm
6110EU		Vari-Lens™ microscope drape for Zeiss™, featuring rotational angled lens	10	122 x 203 cm
6110VL2		Microscope Drape, Zeiss™ with Vari-Lens2	10	122 x 203 cm




## MICROSCOPE DRAPES FOR ZEISS™ MD MICROSCOPES (65 MM LENS)

Zeiss

Part. no	Sketch	Description	Qty/cs	Size
3717UK		Head cover drape for Zeiss™ MD microscope with lens	40	61 cm circular

## MICROSCOPE DRAPES FOR ZEISS™ OPMI MICROSCOPES (48 MM LENS)

Zeiss

Part. no	Sketch	Description	Qty/cs	Size
4941UKCL		Microscope drape for Zeiss™ Opmi™, featuring Clearlens™	20	104x163 cm
4940UK 4940CL		Microscope drape for Zeiss™ Opmi™ Microscope drape for Zeiss™ Opmi™, featuring Clearlens™	10	104x267 cm
2938CL		Microscope drape for Zeiss™ Opmi™, with 2 extended observing oculars	10	104x305 cm

## RECOMMENDED MICROSCOPE DRAPES FOR LEICA™ MICROSCOPES



Leica Microscope Body	Suggested Drapes			
<b>Ophthalmology</b>	Flat Lens	ClearLens	Armatec Angled Lens	VariLens2
M620	2952UK	2952CLEU	8033652EU	6150VL2
M690, M691, M840, M841, M844, M820, M822, Proveo 8	3548EU			6150VL2
<b>ENT</b>	Flat Lens	ClearLens	Armatec Angled Lens	VariLens2
M300, M320, M651, M655, M400 E, M690	2952UK	2952CLEU		6150VL2
MS-2, MS-3	2952UK	2952CLEU		6150VL2
M525 (F20 - floorstand)			8033652EU	6150VL2
M525 MC1			8033651EU	6140VL2
M525 CT20 Ceilingmount			8033652EU	6150VL2
M530 OH6 / OHX			8033651EU	6140VL2
M695	5952UK	5952UKCL		6140VL2
<b>Neurosurgery</b>	Flat Lens	ClearLens	Armatec Angled Lens	VariLens2
M520/M525 OH4			8033651EU	6140VL2
M525 with FL400			8033621EU	
M525 (F50 - Floorstand)			8033651EU	6140VL2
M525 (C50 - Ceilingmount)			8033652EU	6150VL2
M720 OH5 (also when using fluorescence)			8033654EU	
M530 OH6, OHX, Arveo, Provideo			8033651EU	6140VL2
<b>Plastics</b>	Flat Lens	ClearLens	Armatec Angled Lens	VariLens2
M680 without ringmount	3598UK			
M680 with ringmount			8033621EU	
M520/M525 (F20, F40, F50 - floorstand)	2952UK	2952CLEU	8033651EU	6150VL2
M530 OH6 / OHX			8033651EU	6140VL2

For more information about the drape sizes and packed quantities see page 6 and 7.

## LEICA DRAPE SIZES AND PACKED QUANTITIES



Part. no	Sketch	Description	Qty/cs	Size
2952UK 2952CLEU		Microscope drape for Leica™ M690/M655 Microscope drape for Leica™ M690/M655, featuring Clearlens™	10	104x163 cm
8033652EU		Armatec™ microscope drape for Leica™ OH3, OH4, OHS, OH1, MS-2, MS-3, M520/M525, with drawcord, featuring angled lens with click closure	10	117x203 cm
5986UK 5986CLEU		Microscope drape for Leica™ Microscope drape for Leica™, featuring Clearlens™	10	104x267 cm
6150VL2		Microscope Drape, Leica™ with Vari-Lens2	10	117x203 cm
3548EU		Microscope Drape for Leica, with foam lens cover	10	117x267 cm
5952UK 5952UKCL		Microscope drape for Leica™ M690/MS-1 Microscope drape for Leica™ M690/MS-1, featuring Clearlens™	10	137x267 cm

Continues on next page

## LEICA DRAPE SIZES AND PACKED QUANTITIES



Part. no	Sketch	Description	Qty/ cs	Size
3598UK		Microscope drape for Leica™ M680, lens diameter 10 cm	10	117x305 cm
8033651EU 8033650EU		Armatec™ microscope drape for Leica™ OH3, OH4, OHS, OH1, MS-2, MS-3, M520/M525, featuring angled lens with click closure Armatec™ microscope drape for Leica™ OH3, OHS, OH1, MS-2, MS-3, M520/M525, featuring angled lens with click closure	10	117x305 cm 137x381 cm
8033654EU		Armatec™ microscope drape for Leica™ M720 featuring angled lens	10	117x305 cm
8033621EU		Armatec™ microscope drape for Leica™ M680 with ringmount, and M525 with FL400, featuring angled lens	10	117x305 cm
6130VL2		Microscope Drape, Leica™ with Vari-Lens2	10	117x305 cm
6140EU		Vari-Lens™ microscope drape for Leica™, featuring rotational angled lens	10	132x381 cm
6140VL2		Microscope Drape, Leica™ with Vari-Lens2	10	132x381 cm
6135EU		Vari-Lens™ microscope drape for Leica™, featuring rotational angled lens	10	132x305 cm
3652UK 3652CL		Microscope drape for Leica™ (OHS/OH-1 or Wild™ Space Pointer) Microscope drape for Leica™ (OHS/OH-1 or Wild™ Space Pointer), featuring Clearlens™	10	137x381 cm

## MICROSCOPE DRAPES FOR LEICA™ MICROSCOPES - II



Part. no	Sketch	Description	Qty/ cs	Size
890612		Head cover drape for Leica™ microscope with lens	15	56 cm circular

## LASER DRAPES

Part. no	Sketch	Description	Qty/ cs	Size
2903UK		Microscope laser drape, with attached laser arm drape	10	104x305 cm

## HANDLE COVERS

Part. no	Sketch	Description	Qty/ cs	Size
3260UK		Universal handle covers	20	15x36 cm

# ECOLAB AIMS TO REDUCE HAIs & IMPROVE PATIENT AND STAFF SAFETY

We provide innovative solutions to help reduce the risks of infectious disease transmission; to help protect healthcare workers in every type of clinical environment; and to help improve clinical outcomes by maximizing the protection available for both patients and medical professionals during procedures.




## ECOLAB® OPERATING ROOM PROGRAM

### About Ecolab Healthcare

Ecolab Healthcare is driven to help health systems and hospitals realise clinical, operational and financial value through repeatable and measurable workflows. Our standardised processes and digital dashboards provide actionable insights and opportunities for corrective actions that help reduce the costs and inefficiencies of infections, while improving margins and keeping patients and staff safe.

### Together, we can focus on solving your daily challenges

For further information and to discuss a trial, please contact your local Account Executive or Clinical Implementation Specialist.



## ONEMED

Markkinointi: OneMed Oy  
PL 10 (Metsäläntie 20), 00321 Helsinki  
p. 020 786 68  
asiakaspalvelu@onemed.com  
www.onemed.fi

**ECOLAB  
HEALTHCARE**  
EUROPE  
Richtistrasse 7  
8304 Wallisellen  
Switzerland



**MICROTEK MEDICAL, INC.**  
1 Ecolab Place  
St. Paul, MN 55102 U.S.A.  
www.microtekmed.com



**MICROTEK MEDICAL B.V.**  
(a division of Ecolab)  
Hekkehorst 24  
NL-7207 BN Zutphen  
The Netherlands



**MICROTEK MEDICAL  
MALTA LTD.**  
Sorbonne Centre F20,  
Mosta Technopark,  
MST3000, Malta



Microtek Medical B.V.  
Hekkehorst 24  
7207 BN Zutphen  
The Netherlands



www.ecolab.com